

Standard Single Flowmeter Oxygen 0-30Lpm

Product Code: **F1601-30**



Specifications:

OXYGEN
Brass Chrome Plated Body
NOT MRI Conditional
9/16" UNF Fitting
British Standard 5682 Probe
4 bar nominal pressure
BS Colour Fine adjustment Control
Tubing Nipple
0-30Lpm (O₂)
7 Year Manufacturers Warranty

Weights/Dimensions:

Weight: 0.30kg
Height: 204mm
Width: 32mm
Depth: 149mm

Accessories:

SP46 - Tubing Nipple female 9/16 UNF (White for O₂)
H110 - Humidifier (Bubble Type) Autoclavable 250ml Bottle
IS1174 - Oxygen Connecting Tube (1.8m)
SKF6-02 - 5 Year Service Pack

Designed around the needs of Healthcare and Emergency Service Providers, Oxylitre's range of Standard Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Standard Medical Pipeline Flowmeter Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a **7 year** warranty