

## Regulator 3/8bsp Outlet Male Air (MA-4)

Product Code: **R16105A**



### Specifications:

Medical Air  
Brass Chrome Plated Body  
NOT MRI Conditional  
Pin-Index Cylinder Fitting  
138 bar maximum inlet pressure  
Pre-Set 4 bar output pressure  
High Pressure Cylinder Content Gauge  
3/8" Nut & Liner Connection  
**7 Year** Manufacturers Warranty

### Weights/Dimensions:

Weight: 1.25 kg  
Height: 189mm  
Width: 96mm  
Depth: 177mm

### Accessories:

OX010 - Pin-Index "BODOK" Seal (each)  
SK60 - Service Kit for R16 Regulator Range (Bullnose & Pin-Index)

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a **7 year** warranty.