

Standard Single Flowmeter Oxygen 0-15Lpm

Product Code: **F1601**



Specifications:

OXYGEN

Brass Chrome Plated Body

NOT MRI Conditional

9/16" UNF Fitting

British Standard 5682 Probe

4 bar nominal pressure

BS Colour Fine adjustment Control

Tubing Nipple

0-15Lpm (O₂)

7 Year Manufacturers Warranty

Weights/Dimensions:

Weight: 0.30kg

Height: 204mm

Width: 32mm

Depth: 149mm

Accessories:

SP46 - Tubing Nipple female 9/16 UNF (White for O₂)

H110 - Humidifier (Bubble Type) Autoclavable 250ml Bottle

IS1174 - Oxygen Connecting Tube (1.8m)

SKF6-02 - 5 Year Service Pack

Designed around the needs of Healthcare and Emergency Service Providers, Oxylitre's range of Standard Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Standard Medical Pipeline Flowmeter Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a **7 year** warranty