



Standard Single Flowmeter Oxygen 0-15Lpm complete with Humidifier Bottle

Product Code: **F1601H**



Specifications:

OXYGEN
Brass Chrome Plated Body
NOT MRI Conditional
9/16 UNF Fitting
British Standard 5682 Probe
4 bar nominal pressure
BS Colour Fine adjustment Control
Humidifier Bottle complete with Tubing Nipple
0-15Lpm (O₂)
Humidifier Bubble Type:
250ml Jar Autoclavable to 121C
Min and Max Fill Positions
9/16 UNF DISS Fitting (Compatible with Oxygen and Air)
Optimum use between 2lpm and 6lpm
4 bar nominal pressure
Sintered Bronze filtration
7 Year Manufacturers Warranty

Weights/Dimensions:

Weight: 0.38 kg
Height: 340mm
Width: 75mm
Depth: 145mm

Accessories:

H110 - Humidifier (Bubble Type) Autoclavable 250ml Bottle
IS1174 - Oxygen Connecting Tube (1.8m)

Designed around the needs of Healthcare and Emergency Service Providers, Oxylitre's range of Standard Medical Pipeline Flowmeters are the Flowmeter of choice for delivering flexible solutions. Our Standard Medical Pipeline Flowmeter Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Medical Pipeline Flowmeters offer high performance with durability engineered in to the robust design.

Oxylitre's range of Standard Medical Pipeline Flowmeters are dependable, extremely easy to maintain with a long service life and a **7 year** warranty