

# Oxylitre

## Compact Regulator & Brass Flowmeter Oxygen 0-15lpm

Product Code: **R1400**



### Specifications:

OXYGEN

Stainless Steel Body

Bull Nose Cylinder Fitting (UK Size F & G Cylinders)

138 bar maximum inlet pressure

Pre-Set 4 nominal bar output pressure

High Pressure Luminous Cylinder Content Gauge

Flowmeter Brass Chrome Plated:

3/8 BSP Nut & Liner Connection

Tubing Nipple

9/16 UNF Male Outlet

0-15Lpm (O<sub>2</sub>)

**7 Year** Manufacturers Warranty

### Weights/Dimensions:

Weight: 0.84 kg

Height: 216mm

Width: 148mm

Depth: 62mm

### Accessories:

SP46 - Tubing Nipple female 9/16 UNF (White for O<sub>2</sub>)

IS1174 - Oxygen Connecting Tube (1.8m)

BS110 - Bullnose O' Ring Standard Type (each)

SK40 - 5 Year Service Pack

Oxylitre' range of Compact Medical Regulators, have been designed with the end user in mind. Our 138 bar rated Compact Range of regulators are simple to use, reliable and extremely functional, Oxylitre' 138 bar Compact Range Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Compact Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.