



Compact Select-A-Flow Regulator 0-15lpm & Self Sealing Valve Air (MA-4) Bullnose

Product Code: **R1425**



Specifications:

- Medical Air 4 bar (MA-4)
- Stainless Steel Body
- Bullnose Cylinder Fitting (Compatible with UK Size F & G Medical Air Cylinders)
- 138 bar maximum inlet pressure
- Pre-Set 4 bar output pressure
- High Pressure Luminous Cylinder Content Gauge
- Complete with Flowmeter:
 - Select-A-Flow Flowmeter (11 Settings)(Other flow rates available)
 - Tubing Nipple
 - 0-15Lpm (MA-4)
- BS Gas Index Coloured Self Sealing Valve
- 7 Year** Manufacturers Warranty

Weights/Dimensions:

- Weight: 0.90 kg
- Height: 151mm
- Width: 94mm
- Depth: 72mm

Accessories:

- IS1174 - Oxygen Connecting Tube (1.8m)
- BS110 - Bullnose O' Ring Standard Type (each)
- SK40 - 5 Year Service Pack

Oxylitre' range of Compact Medical Regulators, have been designed with the end user in mind. Our 138 bar rated Compact Range of regulators are simple to use, reliable and extremely functional, Oxylitre' 138 bar Compact Range Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Compact Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.