

## Standard Regulator & Brass Flowmeter Oxygen 0-15lpm

Product Code: **R1600**



### Specifications:

OXYGEN

Brass Chrome Plated Body

NOT MRI Conditional

Bullnose Cylinder Fitting

(Compatible with UK Size F & G Oxygen Cylinders)

138 bar maximum inlet pressure

Pre-Set 4 bar output pressure

High Pressure Cylinder Content Gauge

Complete with Flowmeter:

Flowmeter Brass Chrome Plated 9/16 UNF

3/8" Nut & Liner Connection

Tubing Nipple

0-15Lpm (O<sub>2</sub>)

**7 Year Manufacturers Warranty**

### Weights/Dimensions:

Weight: 1.28 kg

Height: 228mm

Width: 183mm

Depth: 76mm

### Accessories:

SP46 - Tubing Nipple female 9/16 UNF (White for O<sub>2</sub>)

IS1174 - Oxygen Connecting Tube (1.8m)

BS110 - Bullnose O' Ring Standard Type (each)

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a **7 year** warranty.