

Regulator 11/16" Nut & Tubing Outlet Carbon Dioxide

Product Code: **R16105-C**



Specifications:

Carbon Dioxide
Brass Chrome Plated Body
NOT MRI Conditional
Pin-Index Cylinder Fitting
50 bar inlet pressure
Pre-Set 4 bar output pressure
Cylinder content Gauge
CO2 11/16" Nut & Tubing Outlet Connection
7 Year Manufacturers Warranty

Weights/Dimensions:

Weight: 1.29 kg
Height: 189mm
Width: 132mm
Depth: 177mm

Accessories:

OX010 - Pin-Index "BODOK" Seal (each)
SK60 - Service Kit for R16 Regulator Range (Bullnose & Pin-Index)

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a **7 year** warranty.