

Regulator High Pressure Medical Air (MA-7)

Product Code: **R1700**



Specifications:

Medical Air 7 Bar (MA-7)
Brass Chrome Plated Body
NOT MRI Conditional
Bullnose Cylinder Fitting
(Compatible with UK Size F & G Oxygen Cylinders)
138 bar maximum inlet pressure
Pre-Set 7 bar output pressure
High Pressure Cylinder Content Gauge
Low Pressure Gauge 0 to 150 PSI
Adjustable Output Pressure
British Standard Gas Indexed Coloured Self Sealing Valve
7 Year Manufacturers Warranty

Weights/Dimensions:

Weight: 1.59 kg
Height: 154mm
Width: 156mm
Depth: 135mm

Accessories:

AH103-3 - Extension Hose Assembly BS Probe to SSV
BS110 - Bullnose O' Ring Standard Type (each)

Oxylitre's range of Standard Medical Regulators, have been designed with the end user in mind. Our Standard Range have a long and exemplary history, simple to use, reliable and extremely functional, our Standard Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Standard Medical Regulators are dependable, easy to maintain with a long service life with a **7 year** warranty.