

## Compact Regulator & Brass Flowmeter Complete With Humidifier Bottle Oxygen 0-15lpm

Product Code: **R2400H**



### Specifications:

OXYGEN

Stainless Steel Body

Bullnose Cylinder Fitting

(Compatible with UK Size F & G Oxygen Cylinders)

200 bar maximum inlet pressure

Pre-Set 4 bar output pressure

High Pressure Luminous Cylinder Content Gauge

Complete with Flowmeter:

Flowmeter Brass Chrome Plated 9/16 UNF

Humidifier Bottle (Autoclavable to 121°C) complete with

Tubing Nipple

0-15Lpm (O<sub>2</sub>)

**7 Year** Manufacturers Warranty

### Weights/Dimensions:

Weight: 0.93 kg

Height: 330mm

Width: 148mm

Depth: 62mm

### Accessories:

H110 - Humidifier (Bubble Type) Autoclavable 250ml Bottle

IS1174 - Oxygen Connecting Tube (1.8m)

BS110 - Bullnose O' Ring Standard Type (each)

SK40 - 5 Year Service Pack

Oxylitre' range of Compact Medical Regulators, have been designed with the end user in mind. Our 200 bar rated Compact Range of regulators are simple to use, reliable and extremely functional, Oxylitre' 200 bar Compact Range Regulators offer high performance with durability engineered in to the robust design.

Oxylitre' range of Compact Medical Regulators are dependable, easy to maintain with a long service life with a 7 year warranty.