

## "elite" On/Off Pipeline Suction Controller Low Suction Direct British Standard Probe

Product Code: **S724**



### Specifications:

- Low Vacuum High Flow (Available in High and Thoracic)
- Lightweight Robust Moulded Construction
- Black Colour Coded Body - Low Suction
- Integral Gauge colour coded - Low Suction
- 0-200mmHg: DUAL Scale: mmHg & kPa (0-26kPa)
- Positive Pressure Relief Valve
- On/Off Flag Style Lever
- Fine adjustment On Off Increase/Decrease Control
- British Standard 5682 Probe (Various Connection Options)
- Encapsulated Hydrophobic/Bacterial Filter Unit
- 7 Year** Manufacturers Warranty
- MRI Compatible

### Weights/Dimensions:

- Weight: 0.52 kg
- Height: 175mm
- Width: 92mm
- Depth: 143mm

### Accessories:

- S750 - Hydrophobic/Bacterial Filters (Box of 30)
- S7500A - 2 Litre Suction Jar Assembly "V" Mounted
- S7500B - 2 Litre Suction Jar Assembly Rail Mounted
- SK17 - 5 Year Service Pack

**Note:** Available in a Twin Low/Low and Twin Low/High Variant

Oxylitre's range of Medical Pipeline Units deliver flexible range of suction solutions. Our Medical Pipeline Suction Range have a long and exemplary history, simple to use, reliable and extremely functional. The Medical Pipeline Suction Range offers high performance with operational reliability engineered in to the robust design.

Oxylitre's range of On/Off Medical Pipeline Suction are dependable, with a minimum maintain required. Our Range of Medical Pipeline Suction Units come with a long service life and a **7 year warranty**.