



"elite" On/Off Pipeline Suction Controller Low Suction Direct British Standard Probe

Product Code: **S724R**



Specifications:

- Low Vacuum High Flow
- Robust Moulded Construction
- Integral Gauge colour coded 0-200mmHg: DUAL Scale: kPa & mmHg
- Positive Pressure Relief Valve
- On/Off Flag Style Lever
- Fine adjustment On Off Increase/Decrease Control
- Encapsulated Bacterial/Hydrophobic Filter Unit
- Complete With Universal Rail Clamp
- 1.5 metre BS Coloured Hose and BS Probe connection
- 7 Year** Manufacturers Warranty
- MRI Compatible

Weights/Dimensions:

- Weight: 0.93 kg
- Height: 175mm
- Width: 92mm
- Depth: 160mm

Accessories:

- S750 - Hydrophobic/Bacterial Filters (Box of 30)
- S7500A - 2 Litre Suction Jar Assembly "V" Mounted
- S7500B - 2 Litre Suction Jar Assembly Rail Mounted
- SK17 - 5 Year Service Pack

Note: Available in a Twin Low/Low and Twin Low/High Variant

Oxylitre's range of Medical Pipeline Units deliver flexible range of suction solutions. Our Medical Pipeline Suction Range have a long and exemplary history, simple to use, reliable and extremely functional. The Medical Pipeline Suction Range offers high performance with operational reliability engineered in to the robust design.

Oxylitre's range of On/Off Medical Pipeline Suction are dependable, with a minimum maintain required. Our Range of Medical Pipeline Suction Units come with a long service life and a **7 year warranty**.