

## "elite" On/Off Pipeline Suction Controller Thoracic Suction Direct British Standard Probe

Product Code: **S725**



### Specifications:

Thoracic Vacuum High Flow (Available in High and Low)  
Lightweight Robust Moulded Construction  
Green Colour Coded Body - Thoracic Suction  
Integral Gauge colour coded - Thoracic Suction  
0-70cmH<sub>2</sub>O: DUAL Scale: cmH<sub>2</sub>O & kPa (0-5kPa)  
Positive Pressure Relief Valve  
On/Off Flag Style Lever  
Fine adjustment On Off Increase/Decrease Control  
British Standard 5682 Probe (Various Connection Options)  
Encapsulated Hydrophobic/Bacterial Filter Unit  
**7 Year** Manufacturers Warranty  
MRI Compatible

### Weights/Dimensions:

Weight: 0.52 kg  
Height: 175mm  
Width: 92mm  
Depth: 143mm

### Accessories:

S740 - Bacterial Suction Filters (Box of 30)  
S7500A - 2 Litre Suction Jar Assembly "V" Mounted  
S7500B - 2 Litre Suction Jar Assembly Rail Mounted  
SK17 - 5 Year Service Pack

**Note:** Available in a Twin High/Thoracic Variant

Oxylitre's range of Medical Pipeline Units deliver flexible range of suction solutions. Our Medical Pipeline Suction Range have a long and exemplary history, simple to use, reliable and extremely functional. The Medical Pipeline Suction Range offers high performance with operational reliability engineered in to the robust design.

Oxylitre's range of On/Off Medical Pipeline Suction are dependable, with a minimum maintain required. Our Range of Medical Pipeline Suction Units come with a long service life and a **7 year warranty**.