

Entonox® Pain Relief Demand Valve (Mk2)



ENTONOX® DEMAND UNIT & REGULATOR

Completed with:

- 1 x Adult Face Mask (disposable type)
- 5 x Mouthpieces
- 1 x Bacterial Filter
- 1 x Pin-Index Cylinder Key

PRU2001
PRU2001L

Hose Length 1.5m
Hose Length 3.0m



ENTONOX® DEMAND UNIT & REGULATOR IN CARRY BAG

Completed with:

- 1 x Adult Face Mask (disposable type)
- 1 x Child Face Mask (disposable type)
- 5 x Mouthpieces
- 1 x Bacterial Filter
- 1 x Pin-Index Cylinder Key

PRU2002
PRU2002L

Hose Length 1.5m
Hose Length 3.0m



ENTONOX® DEMAND UNIT & BS PROBE ONLY

Completed with:

- 1 x Adult Face Mask (disposable type)
- 5 x Mouthpieces
- 1 x Bacterial Filter

Also available for connecting direct to piped ENTONOX® gas systems.

(Longer hose lengths available on request)

PRU2003
PRU2003L

Hose Length 1.5m
Hose Length 3.0m



PRU050

- Carry Case only
- Luminous Stitching
- AEGIS Antimicrobial Technology



PRU01035

- Bacterial Filter & Mouthpiece (50 per box)



PRU2003/W

- Demand Unit Wall Mounting Bracket



PRU0-21

- Demand Unit Waste Gas Outlet



FACE MASKS

	Disposable	Silicone (re-usable)
Adult	011-4	012-4
Child	011-2	012-2

Specifications:

- 7 Year Warranty
- 5 Year Service Intervals
- Exhaust Vent Away From User
- Antimicrobial Resistance Infused
- 360° Rotating Head
- Low Patient Effort
- Uses Standard In-Line Filters
- Available in PVC or Silicone Hose
- Ergonomically Designed Handle
- Robust, Lightweight Body

Oxylitre's innovated PRU demand valve is manufactured to the highest possible quality, designed for a minimum of 10+ year lifecycle. The unit offers the end user full 360° rotation around the hose, also incorporating a sleek handle for a comfortable patient experience. Set within a robust, lightweight, moulded body and complete with a British Standard Stainless Steel probe, our innovative engineers have designed this with the patient in mind. For more information, or to discuss pricing & service ex, please contact your local Oxylitre representative